


3M™ Littmann® Master Classic II Stethoscope



3M Medical
from the 3M Health Care family



3M™ Littmann® Master Classic II Stethoscope



The excellent acoustic sensitivity of a tunable diaphragm is the hallmark of the low-profile Littmann Master Classic II stethoscope. Award-winning design combines high performance with portable convenience in this single-sided stethoscope.

Features

- Tunable diaphragm conveniently alternates between low- and high-frequency sounds without the need to turn over the chestpiece
- Excellent acoustic seal and comfortable fit with patented 3M™ Littmann® Snap Tight Soft-Sealing eartips
- Portable convenience of a single-sided chestpiece
- Comfortably angled, anatomically correct headset
- Patient-friendly, nonchill rim and diaphragm
- High-performance, single-lumen tubing in a range of attractive colors
- Three-year warranty

Ordering Information

Catalog #	Length	Tube Color
2144L	27 in. (69cm)	Black
2143	27 in. (69cm)	Purple
2146	27 in. (69cm)	Burgundy
2147	27 in. (69cm)	Navy Blue
2630	27 in. (69cm)	Caribbean Blue
2632	27 in. (69cm)	Hunter Green
2633	27 in. (69cm)	Ceil Blue
<i>Black Edition — Black-plated chestpiece and eartubes</i>		
2141	27 in. (69cm)	Black
<i>Gold Edition — 23-Karat gold-plated chestpiece and eartubes</i>		
2142G	27 in. (69cm)	Black

The Littmann stethoscopes listed do not contain natural rubber latex or dry natural rubber as components in the product or packaging.

For more information, call the 3M Health Care Helpline at **1-800-228-3957**, or visit us at www.littmann.com



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957

www.littmann.com

3M Canada Inc.

Post Office Box 5757
London, Ontario N6A4T1
Canada
1 800-563-2921

Outside of USA, please contact your 3M Representative



40% Pre-consumer waste paper
10% Post-consumer waste paper

Printed in U.S.A.
Copyright © 3M 2005
All Rights Reserved
70-2008-4251-9